



THE CENTER FOR
GENOMICS AND
BIOINFORMATICS

Emily McKillip and Kris Klueg

DGRC Vector Standard Operating Procedure 1.0:
Processing Clones for Whatman[®] FTA[®] Discs

CGB Technical Report 2006-02

doi:10.2506/cgbtr-200602

McKillip, E. and Klueg, K. 2006. DGRC Vector Standard Operating Procedure 1.0: Processing Clones for Whatman[®] FTA[®] Discs. *CGB Technical Report 2006-02*. The Center for Genomics and Bioinformatics, Indiana University, Bloomington.

**DGRC Vector Standard Operating Procedure 1.0: Processing Clones for
Whatman® FTA® Discs**
(adapted from the Whatman protocol)

Emily McKillip and Kris Klueg

Drosophila Genomics Resource Center, Center for Genomics and Bioinformatics, Indiana
University, Bloomington, Indiana, 47405

Please process clones within a month of receipt for best results:

Items needed:

- 1X sterile TE
- Thawed competent cells on ice -
 - * We recommend DH5alpha or an equivalent strain
 - * Cells should have a transformation efficiency of 1×10^7 cfu/ug or higher
- 1 ml LB media aliquots
- LB (+ appropriate antibiotic) plates (see p. 2 for additional considerations)

To process:

1. Add 50 μ l of 1X sterile TE to the microfuge tube with the clone disc and pipet up and down quickly 2X - REMOVE the "TE" IMMEDIATELY. This is the critical step – **do not leave the TE on the disc for more than one or two seconds (see note below).**

Please note: The TE step is necessary to wash off the chemical that Whatman has placed on the filter – it is this chemical that facilitates the lysis of the bacterial clones and allows the DNA to stick. Failure to do the TE wash step will inhibit transformation. However, it is important to note that if the TE sits on the disc for more than one or two seconds, the DNA will elute off and be lost, resulting in a failed transformation.

2. Place tube on ice
3. Add 50 μ l of competent cells and incubate on ice for 30 minutes. Transformation efficiency will increase if you vortex the competent cells/disc mixture half way through the ice incubation step and once more before the heat shock. Vortexing should be for one second only and return immediately to ice.
4. Heat shock cells/disc for 2 minutes at 37°C
5. Transfer cells only (leave disc behind) to 1 ml of LB media and incubate with shaking at 37°C for one hour

6. Plate 100 μ l and 200 μ l on LB (+antibiotic) plates. **For GOLD and EST clones we recommend plating the entire transformation.**

Ideally you will obtain approximately 1 colony for every microliter of the transformation plated. However, there will be considerable variability among clones.

If you have followed these procedures exactly and did not obtain colonies, please contact Kris Klueg (dgrc@cgb.indiana.edu).

Other important considerations:

Subcloning efficiency competent cells are **not** recommended. If you are purchasing cells commercially, we recommend library efficiency. If you are making your own cells, the transformation protocol on our website typically gives cells that are 1×10^7 cfu/ μ g.

For chloramphenicol resistant clones, do not use XL10 GOLD competent cells. These cells carry a chloramphenicol-resistance gene on the F' episome

Chloramphenicol – the concentration requirement in plates (for transformation) versus media (for large scale growth) differs. Please see our FAQ page on our website for the concentrations

cDNAs from the RE library (REXXXXXX) give low DNA yields when grown on ampicillin. For large scale DNA preps, we recommend you use carbenicillin (an ampicillin analog) at 75 μ g/ml. (Using carbenicillin for the initial transformation is not necessary).

We recommend using carbenicillin for large scale DNA preps of many of the common vectors (pUAST derivatives, pCaSpeR derivatives, and Gateway).

Acknowledgements:

This project was supported by Grant Number 1 P40 RR017093 from the National Center for Research Resources (NCRR), a component of the National Institutes of Health (NIH). Its contents are solely the responsibility of the authors and do not necessarily represent the official views of NCRR or NIH.